



Phillips Respironics EverGo™



Weight: 8.5 lb
Dims: 12" x 6" x 8.5"

1 Battery / with 2nd Battery
Setting 1 = 6 Hr / 12 Hr
Setting 2 = 4 Hr / 8 Hr
Setting 3 = 3 Hr / 6 Hr
Setting 4 = 2 Hr / 4 Hr
Setting 5 = 2 Hr / 4 Hr
Setting 6 = 2 Hr / 4 Hr



Weight: 9.7 lb
Dims: 11.62" x 6" x 12.39"

Inogen One™

Setting 1 = 3.0 Hr
Setting 2 = 3.0 Hr
Setting 3 = 2.5 Hr
Setting 4 = 2.5 Hr
Setting 5 = 2.0 Hr



Weight: 4.9 lb
Dims: 9.5" x 7.5" x 3.125"

Inova Labs International Biophysics™ LifeChoice

Internal / with External Battery
Setting 1 = 2 Hr / 5 Hr
Setting 2 = 2 Hr / 5 Hr
Setting 3 = 2 Hr / 5 Hr



Weight: 6 lb
Dims: 10" x 7" x 4"

Invacare XPO2 - XPO100™

Internal / with External Battery
Setting 1 = 3.5 Hr / 7 Hr
Setting 2 = 2.5 Hr / 5 Hr
Setting 3 = 2.0 Hr / 5 Hr
Setting 4 = 1.5 Hr / 3 Hr
Setting 5 = 1.0 Hr / 2 Hr



Weight: 4.4 lb
Dims: 8.6" x 6.1" x 3.6"

AirSep FreeStyle™

Internal Battery / with Battery Belt
Setting 1 = 3.5 Hr / 10 Hr
Setting 2 = 2.5 Hr / 6 Hr
Setting 3 = 2.0 Hr / 5 Hr



Weight: 9.75 lb
Dims: 5.5" x 7.25" x 16.31"

AirSep LifeStyle™

Setting 1 = 3.25 Hr*
Setting 2 = 3.25 Hr*
Setting 3 = 3.25 Hr*
Setting 4 = 3.25 Hr*
Setting 5 = 3.25 Hr*
*with Power Pack (4 batteries)
Each Battery lasts 50 min



Weight: 9.8 lb
Dims: 11.6" x 7.4" x 4.6"

Delphi Central Air, Evo™

Setting 1 = 3 Hr 25 min
Setting 2 = 2 Hr 45 min
Setting 3 = 2 Hr 15 min
Setting 4 = 1 Hr 40 min
Setting 5 = 1 Hr 15 min



Weight: 7 lb
Dims: 11.61" x 5.98" x 10.75"

Inogen One G2™

Standard / Extended
Setting 1 = 4.0 Hr / 8.0 Hr
Setting 2 = 3.5 Hr / 7.0 Hr
Setting 3 = 2.7 Hr / 5.4 Hr
Setting 4 = 2.2 Hr / 4.4 Hr
Setting 5 = 2.0 Hr / 4.0 Hr



Weight: 19.9 lbs
Dims: 16.5" x 11" x 8"

Invacare Solo™

Pulse / Continuous
Setting .5 = N/A / 4.5 Hr
Setting 1 = 4.5 Hr / 3.5 Hr
Setting 2 = 3.5 Hr / 2.5 Hr
Setting 3 = 3.0 Hr / 1.5 Hr
Setting 4 = 2.5 Hr
Setting 5 = 2.5 Hr



Weight: 17.9 lb
Dims: 19.3" x 12.3" x 7.1"

SeQual Eclipse®

Pulse / Continuous
Setting .5 = NA / 4.4 Hr
Setting 1 = 4.4 Hr / 3.7 Hr
Setting 2 = 3.6 Hr / 2.4 Hr
Setting 3 = 3.0 Hr / 1.3 Hr
Setting 4 = 2.6 Hr
Setting 5 = 2.3 Hr
Setting 6 = 2.1 Hr



Weight: 19 lb
Dims: 15" x 11" x 8"

DeVilbiss Healthcare iGo™

Pulse / Continuous
Setting 1 = 5.4 Hr / 4.0 Hr
Setting 2 = 4.7 Hr / 2.4 Hr
Setting 3 = 4.0 Hr / 1.6 Hr
Setting 4 = 3.5 Hr
Setting 5 = 3.2 Hr
Setting 6 = 3.0 Hr



Weight: 14.85 lb
Dims: 8" x 8" x 12"

OxLife Independence™

Pulse / Continuous
Setting 1 = 2.8 Hr / 2.4 Hr
Extended Battery = 5.6 Hr / 4.8 Hr
Setting 2 = 2.8 Hr / 1.4 Hr
Extended Battery = 5.6 Hr / 2.8 Hr
Setting 3 = 1.4 Hr / 53 min
Extended Battery = 2.8 Hr / 1.8 Hr
Setting 4 = 1.4 Hr
Extended Battery = 2.8 Hr
Setting 5 = 53 min
Extended Battery = 1.8 Hr
Setting 6 = 53 min
Extended Battery = 1.8 Hr

Battery Duration Chart Estimated Times Based on NEW Batteries Hr = hours / min = minutes / lb = pounds / Dims = dimensions

POC Short Term Rentals, Cruise Ship Delivery, Hotels, Airport Meets