E OXYGE



ResMed MobiTM Weight: 5.5 **l**bs Dims: 9.6" x 7.95" x 4.49



EasyPulse POC3 Weight: 4.9 lbs Dims: 8.5" x 4.5" x 6.5



Precision Medical Active5 Weight: 5 lbs Dims: 8.4"x3.2"x8.5"



Platinum Mobile Weight: 4.98 lbs Dims: 9.45" x 7.5" x 3.88"

SUPPORT

NEW POC'S WITH FAA APPROVED LABEL



GCE Healthcare Zen-O Weight: 10.3 lbs. Dims: 8.3" x 6.6" x 12.3"

The manufacturer of this POC has determined this device conforms to all applicable FAA requirements for POC carriage and use on board aircraft.

EXAMPLE FAA APPROVED LABEL







Drive DeVilbiss iGo2 Weight: 5 lbs Dims: 8.4"x3.5"x8.6"



Weight: 6.6 lbs Dims: 10.1" x 6.5" x 4.5



Weight: 4.7 lbs Dims: 7.19" x 3.26" x 8.15" (Also the OxyGo Next)



Inogen G4TM Weight: 2.8 lbs. Dims: 7.2" x 5.91" x 2 (Also the OxyGo Fit)



Inogen G3TM Weight: 4.8 **I**bs Dims: 8.75" x 3.0" x 7.25 (Also the OxyGo)



Philips Respironics SimplyGo Weight: 10 lbs. Dims: 11.5"x 10"x 6"



Philips Respironics SimplyGo Mini Weight: 5 lbs. Dims: 10.2" x 8.3" x 2.6'



O2 Concepts Oxlilfe Freedom Weight: 5.5 lbs Dims: 9.4"x9"x3.4



Belluscura XPLO2R Weight: 3.25 lbs Dims: 4cell 7.3"x2.9"x7" Dims: 8cell 7.3"x2.9"x7.6"



Weight: 5 lbs Dims: 8.4"x3.5"x8.6"



куthm Healthcare Weight: 4.37 lbs Dims: 8.7"x3.35"x6.3"



Invacare XPO2TM Weight: 6 lbs Dims: 10" x 7" x 4"



Weight: 11.5 lbs ns: 14.63" x 10.5" x



Weight: 17.9 lbs Dims: 19.3" x 12.3" x 7.1



AirSep Focus TM Weight: 1.75 lbs Dims: 6.2" x 4.6" x 2.4"





Weight: 6.2 lbs Dims: 10.7" x 6.6" x 4.4"



Inova Labs Activ Ox⁷ Weight: 4.83 lbs Dims: 9.05" x 7.875" x 4.38"



Weight: 19 lbs Dims: 15" x 11" x 8'





TALS & SALES

00-346-3556

PORTABLE OXYGEN CONCENTRATORS

AS OF LATEST FAA RULING

THIS DOCUMENT REVISED 2022

AG OF LATEOT FAA HOLING							
FAA APPROVED BATTERY DURATION Setting	0.5	1.0	2.0	3.0	4.0	5.0	6.0
Caire AirSep Focus			1.5 hrs				
Focus w/ Battery AirBelt			4.0 hrs	/			
Caire AirSep Freestyle		6.5 hrs	3.5 hrs	3.0 hrs		(
FreeStyle Battery AirBelt		6.5 hrs	3.5 hrs	3.0 hrs	/	7	
Caire AirSep FreeStyle 5		4.25 hrs	2.5 hrs	2.0 hrs	1.75 hrs	1.25 hrs	
FreeStyle 5 Battery AirBelt		4.0 hrs	2.5 hrs	2.0 hrs	1.75 hrs	1.5 hrs	
Belluscura XPLO2R (8 cell battery)		4hr 44min	4hr 20min	3hr 16min	2hr 20min		
Caire FreeStyle Comfort (Single Battery)		8 hrs	4 hrs	3 hrs	2.25 hrs	2 hrs	
Caire FreeStyle Comfort (Double Battery)		16 hrs	8 hrs	6 hrs	4.5 hrs	4 hrs	
Inogen One G3 8 Cell (Also the OxyGo)		5.0 hrs	4.0hrs	3.0 hrs	2.0 hrs	1.5 hrs	
G3 16 Cell Extended Battery		10.0 hrs	9.0 hrs	6.0 hrs	4.0 hrs	3.0 hrs	
Inogen One G4 Single Battery		2hr 40min	2hr 15min	1hr 20min			
G4 Double Battery (Also the OxyGo Fit)	/	5.0 hrs	4.5 hrs	2.5 hrs	/		
Inogen One G5 Single Battery		6 hrs	4 hrs 45 min	3 hrs	2 hrs	1 hr 30 min	1 hr 15 min
G5 Double Battery (Also the OxyGo Next)	<u> </u>	13 hrs	10 hrs	6 hrs 45 min	4 hrs 30 min	3 hrs 15 min	2 hrs 30 min
Inova Labs LifeChoice Activox		12.0 hrs	6.0 hrs	5.0 hrs			
Activox External Battery +add	-	3.0 hrs	3.0 hrs	2.0 hrs			
Inova Labs Activox Pro		12.5 hrs	7.25 hrs	3.25 hrs :			
Activox External Battery +add	<u> </u>	3.0 hrs	3.0 hrs	3.0 hrs			
Inova Labs Activox 4 L		10.25 hrs	8.25 hrs	5.0 hrs	4.0 hrs		
Activox External Battery +add		3.0 hrs	3.0 hrs	3.0 hrs	3.0 hrs		
Invacare Platinum Mobile		5.08 hrs	3.5 hrs	2.3 hrs	1.75 hrs		
Invacare Platinum Mobile with 2nd Battery		10.17 hrs	7.0 hrs	4.67 hrs	3.5 hrs		
Invacare XPO2		3.67 hrs	2.67 hrs	2.0 hrs	1.67 hrs	1.3 hrs	
XPO2 Ext. Battery +add		3.5 hrs	2.4 hrs	2.0 hrs	1.5 hrs	1.0 hrs	
O2 Concepts Oxlife Freedom 8 Cell		5 hrs	4 hrs	3 hrs	2 hrs	1.75 hrs	
O2 Concepts Oxlife Freedom 12 Cell		10 hrs	8 hrs	6 hrs	4 hrs	2.5 hrs	
Lifestyle P2 Rythm Healthcare		5 hr	4 hr	3 hr	2hr 14min	1hr 50min	
Precision Medical Active5		6.5 hrs	4.3 hrs	2.7 hrs	2 hr	1.5 hrs	
Precision Medical EasyPulse PM4150		4.7 hrs	3.4 hrs	2.7 hrs	2.0 hrs	1.6 hrs	
EasyPulse Ext. Battery +add		4.7 hrs	3.4 hrs	2.7 hrs	2.0 hrs	1.6 hrs	
Precision Medical EasyPulse PM4130		5.5 hrs	3.8 hrs	2.8 hrs	/		
EasyPulse Ext. Battery +add		5.5 hrs	3.8 hrs	2.8 hrs		\	
Philips Respironics SimplyGo Mini	T	6.0 hrs	4.5 hrs	3.0 hrs	2.25 hrs	2.0 hrs	
SimplyGo Mini Double Battery		12.0 hrs	9.0 hrs	6.0 hrs	4.5 hrs	4.0 hrs	
ResMed Mobi Internal		8.6 hrs	6.0 hrs	3.8 hrs	3.2 hrs	/	
ResMed Mobi Ext Battery +add		4.5 hrs	3.0 hrs	1.6 hrs	1.6 hrs	\	
DeVilbiss Healthcare iGo2	<u> </u>	5 hrs	3.5 hrs	2.5 hrs	1.75 hrs	1.5 hrs	
DeVilbiss Healthcare iGo	[5.4 hrs	4.7 hrs	4.0 hrs	3.5 hrs	3.2 hrs	3.0 hrs
Continuous Flow		4.0 hrs	2.4 hrs	1.6 hrs	/ 	\	
Precision Medical EasyPulse TOC]	3.3 hrs	3.0 hrs	2.8 hrs	2.5 hrs	2.3 hrs	
Continuous Flow 1/8th=5.0 hrs 1/4th=4.3 hrs	3.3 hrs	1.8 hrs	0.8 hrs			-	
GCE Healthcare Zen-O]	4.0 hrs	4.0 hrs	3.0 hrs	2.25 hrs	2.0 hrs	1.75 hrs
Zen-O w/ 2nd Battery		8.0 hrs	8.0 hrs	6.0 hrs	4.5 hrs	4.0 hrs	3.5 hrs
Continuous Flow Battery/with 2nd Battery	3.0/6.0hrs	1.75/3.5hrs	0.75/1.5hrs	/-			
Invacare Solo2		4.5 hrs	3.5 hrs	3.0 hrs	2.5 hrs	2.5 hrs	
Continuous Flow		3.5 hrs	2.5 hrs	1.5 hrs	/		
Oxlife Independence	3hr 8min	3.0 hrs	2hr 53min	2.0 hrs	1.5 hrs	1hr 8min	53 min
w/ 2nd Battery	6hr 15min	6.0 hrs	5hr 45min	4.0 hrs	3.0 hrs		1hr 45min
Continuous Flow Battery/with 2nd Battery	2hr53min/5hr45min	2hr38min/5hr15min	1hr15min/2hr30min	45min/1hr30min			l
Philips Respironics SimplyGo		3.4 hrs	3.0 hrs	2.2 hrs	1.7 hrs	1.3 hrs	1.3 hrs
Continuous Flow	2.9 hrs	1.6 hrs	0.9 hrs				
SeQual Eclipse		4.4 hrs	3.6 hrs	3.0 hrs	2.6 hrs	2.3 hrs	2.1 hrs
Continuous Flow	4.4 hrs	3.7 hrs	2.0 hrs	1.3 hrs			

SeQual Eclipse 5

Pulse Delivery
Setting .5 NA
Setting 1.0 5.4 hrs
Setting 2.0 5.2 hrs
Setting 3.0 4.9 hrs
Setting 4.0 4.0 hrs
Setting 5.0 3.7 hrs
Setting 6.0 3.5 hrs

Continuous Flow
Setting 1.0 3.7 hrs
Setting 1.0 3.7 hrs
Setting 2.0 2.0 hrs
Setting 3.0 1.3 hrs
Setting 5.0 3.7 hrs
Pulse 128 mL 2.5 hrs
Setting 6.0 3.5 hrs
Pulse 192 mL 1.7 hrs



All times based on NEW batteries and are estimates as per each manufacturer. Each Airline has its own rules about flying with portable oxygen concentrators. Battery run times are shorter at higher settings.

1-800-346-3556

RENTALS & SALES

24HR SUPPORT

INFO@AEROMEDIC.COM